



Hand Held Calibration / Diagnostic Tool

General features



- Allows easy access to calibration, diagnostic and adjustment menus built into PG Trionic control modules.
- Simplifies setup and testing
- Two line 16 character backlit display
- LEFT / RIGHT buttons select menus and settings to be adjusted
- UP / DOWN buttons adjust settings
- ENTER / ESC move from menu to sub-menu and back
- Complete with 4' replaceable cable

MENU: DIAGNOSTICS	When EZcal connected to card, first menu automatically displayed
DIAGNOSTICS: TILT 0.3', 1.2'	Press ENTER to select the menu; first item in that menu displayed
DIAGNOSTICS: HEIGHT 15%	Press RIGHT to see next item in selected menu
MENU: DIAGNOSTICS	Press ESC to escape from selected menu
MENU: ADJUSTMENTS	Press RIGHT to select a different menu
ADJUSTMENTS: XTILT TRIP 5.0'	Press ENTER to select the new menu; first item in new menu displayed
ADJUSTMENTS: XTILT TRIP 4.2'	Press UP (or DOWN) to adjust the displayed item
ADJUSTMENTS: YTILT TRIP 3.0'	Press RIGHT to select another item